

# HAFFMANS IAM MEDIDOR DE AR INPACK 2000

PRODUCT LEAFLET

## GENERAL PRODUCT INFORMATION

In the beer and beverage industries, the content of dissolved oxygen ( $O_2$ ) is a determining factor in the quality, taste and flavor stability of beer and carbonated beverages, as it causes a rapid decline of the flavor stability and shelf life. By measuring the air volume in the container, the  $O_2$  content can be determined. Breweries as well as soft drink manufacturers continuously control and measure the concentration of air after filling.

The Inpack 2000 Air Meter, type IAM, is used to determine the air content in the headspace volume as well as the total air content in the beverage, in combination with the Inpack 2000  $CO_2$  Calculator or the Inpack 2000  $CO_2$  Meter. The principle is based on  $CO_2$  absorption by caustic.

The pressure of the container is released and the gas fed to the IAM. The  $CO_2$  is completely absorbed by the caustic and the remaining gas is collected in the measuring burette. Its volume (ml) is shown on the measuring burette.

The device is available in two executions:


- Type IAM-20 has a measuring range of 0-20 ml, with higher accuracy in the low range and is typically used for carbonated beverages with low air content such as beer
- Type IAM-36 measures the range of 0-36 ml and is typically used for carbonated beverages with a higher air content, such as soft drinks

## BENEFITS

- Perfect quality management
- Cost saving
  - maintenance free

## APPLICATIONS

- Laboratory, in combination with Pentair Haffmans' Inpack 2000  $CO_2$  Calculator and Inpack 2000  $CO_2$  Meter



# HAFFMANS IAM MEDIDOR DE AR INPACK 2000

PRODUCT LEAFLET

## TECHNICAL DATA

### IAM-20

**Measuring range**

0.0 - 20.0 ml

**Subdivision**

Range 1	0.0 - 0.50 ml: 0.05 ml
Range 2	0.6 - 6.50 ml: 0.10 ml
Range 3	7.0 - 20.00 ml: 0.50 ml

**Accuracy**

± 5 % of measured value  
with a min. of 0.025 ml

**Caustic content**

Approx. 160 ml

**Weight**

Ca. 1.4 kg

### IAM-36

**Measuring range**

0.0 - 36.0 ml

**Subdivision**

Range 1	0.0 - 6.0 ml: 0.10 ml
Range 2	7.0 - 36.0 ml: 0.50 ml

**Accuracy**

± 5 % of measured value  
with a min. of 0.05 ml

**Caustic content**

Approx. 175 ml

**Weight**

Ca. 1.3 kg

**Dimensions**

250 x 170 x 595 (L x W x H mm)

**Caustic concentration**

30 % g/l NaOH with approx. 1 % methanol  
30 % g/l KOH with approx. 1 % methanol


## SCOPE OF SUPPLY

- Inpack 2000 Air Meter\*
  - Type IAM-20 or
  - Type IAM-36
- Measuring burette and leveling vessel with clamp
- Set of tubing
- Instruction manual

\* When ordering, please state which instrument you require.

## OPTIONS

- Certificate of measurement

**HAFFMANS BV**

P.O. BOX 3150 NL-5902 RD VENLO, NETHERLANDS WWW.HAFFMANS.NL

All Pentair trademarks and logos are owned by Pentair, Inc. All other brand or product names are trademarks or registered marks of their respective owners. Because we are continuously improving our products and services, Pentair reserves the right to change specifications without prior notice.

Pentair is an equal opportunity employer.

haffmans\_inpack\_2000\_air\_meter\_en\_1204 © 2012 Pentair, Inc. All Rights Reserved.

